

1

HAND BOOK & INSTRUCTION MANUAL OF Narrow Band UVB Phototherapy Device, 311nm with Magnetic Interchangeable Aperture Plates For Targeted Spot therapy (Hand Held Home Care Device)


Device Includes:

- original Philips NBUBV,311nm Lamp/CE Certified/5-star rated (Made in Poland)
- Magnetic comb attachment
- Interchangeable Magnetic Aperture Plates for targeted spot therapy

A Unique & Innovative Smart Device Helps To Treat Vitiligo/Leucoderma/Safed Daag/ White patches/Fulveri, Psoriasis, Eczema, Atopic dermatitis, Mycosis fungoides, Pruritus, Photodermatoses, Hyperbilirubin and many skin related ailments and complications.



Interchangeable Aperture Plates for Radiation Control (For Various sizes of skin patches)




VITILIGO

BEFORE AFTER



PHOTOMETRIC RESPONSE OF LAMP (311NM)



PSORIASIS

BEFORE AFTER





Vision Tech Electronics was established during the year 2000 in Ahmedabad (Gujarat, India), and attained Trade Mark (Under Trade Registry, Class:9, 10 & 11, Govt of India) entitled: **VISION TECH®**

We have wide experience in manufacturing , more than 20 years in highly specialized Lighting and UV related applications and apparatus. We are very delighted to present Handheld Narrow Band UVB, 311nm along with interchangeable magnetic aperture plates to get desired compatible various patterns of radiation for various sizes of skin patches. hence targeted treatment therapy can be achieved, also separate magnetic comb attachment is provided for the patches on scalp or sensitive skin.

We are using original Philips NBUVB , 311nm , CE Certified Lamp, 5-star rated (Made in Poland). NBUVB device has already helped thousands of the patients across the world in treatment of Vitiligo / Leucoderma / Samedag / White patches / Fulver, Psoriasis, Eczema, Atopic Dermatitis, Mycosis Fungoides, Pruritus, Photodermatoses, Hyperbilirubin and skin related ailments and complications.

Although designed for consumer at-home therapy, also is used by many professional Dermatologists, Clinics and Hospitals.

The precise and unique spectral energy distribution with an emission peak at 311nm, almost nullifies the potential side effects like redness and itching. Also, a shorter period of exposure time is needed.

- Treatment can be carried out alternate day for each area of skin to be treated and in dark room for effective radiation exposure.

- While using the device, simply place it against the affected area (Even you may touch the device). Positive results have been seen in as little as 15 days. however results do vary between individuals and severity of condition. We recommend to consult your Dermatologist for full treatment before use of the device

Our scientific geometrical physical design allows to get 2.5 times more radiation rate as compare to other conventional Phototherapy devices. In our Phototherapy device, mirolite curvature and cavity ensures that Narrowband UVB, 311nm can be safely focused and effectively over the skin to be treated without any deviation.

Interchangeable Aperture Plates:

These have been provided in various sizes to create various radiation patterns, compatible to various patch sizes of the skin.

Magnetic Comb:

We should attach comb part only when radiation exposure is meant for the patches on scalp or sensitive area of skin, otherwise comb part should be avoided. Direct usage of the device (without comb) reduces exposure time and gives the fastest results.

Rare Earth Magnets (Neodymium):

Attachment mechanism between device & Aperture plate + Comb is accomplished by rare earth magnets (Neodymium) which facilitate easy coupling of aperture plates or comb part with the surface of main body of NBUVB device. This magnetic mechanism keeps the stable coupling & keeps away mechanical wear & tear.

These above described features allow the NBUVB Phototherapy sessions in disciplined manner and that systematic exposure helps in fast and targeted recovery from the skin ailments with out any significant side effects and this controlled mechanism avoids unnecessary skin tanning.


Control Gears:

Control gear of Device ensures very long life span of the Lamp with minimal depreciation of exposed energy and has been tested most stringently to run the Lamp very safely. Hence our careful design of NBUVB Phototherapy device helps to give very fast results.

PRODUCT FEATURES


Aesthetic / General Features:

- ➔ Device Body: high quality medical grade virgin ABS material
- ➔ Body colour: Coral White
- ➔ Aperture plates: made of Acrylic (front: Blue, Rear: Mirolite)
- ➔ Comb: made of poly carbonate
- ➔ Rugged & elegant packaging, helps device to keep safe & intact
- ➔ Dimensions (LxBxH): 450mmx60mmx48mm, Weight:450gm



Optical Parameters:

- ➔ Uses original Philips NB-UVB, 311nm CE Certified, 5-Star rated Lamp (Made in Poland)
- ➔ High quality curved mirolite surface beneath the lamp focuses the UV light energy at optimal level
- ➔ Treatment area: 5 cm x 15 cm (2" x 6")
- ➔ Green indicator signifies healthy lamp
- ➔ Blinking Red indicator along with buzzer signifies unhealthy lamp.



Electrical Parameter:

- ➔ Premium quality of industrial Electronic Ballast has been used, which ensures enhanced life span of lamp
- ➔ Mains Supply Voltage: 230Vac
- ➔ Length of 3-pin mains cord: 7ft
- ➔ SPD>6KV, PF>0.90, RFI/EMI Suppression
- ➔ ON/OFF Rocker switch is provided
- ➔ Lamp:PL-S, 9W,Voltage:60V, Current:0.170A, Wavelength:311nm
- ➔ Cape Base: 2-pin, G23, Colour designation: NBUVB

Includes:

High quality Eye Protecting Goggle in separate box & micro fiber have been provided absolutely free and also Instruction manual/ hand book to give systematic & proper guidelines to treat the ailments along with NBUVB Phototherapy device, also highlights the understanding of Vitiligo / Psoriasis.


How to use NBUVB Phototherapy Device:

- ➔ Plug in three pin mains cord into three pin socket, switch on the switch
- ➔ Switch on Rocker switch, which is provided on device
- ➔ Simply keep device against affected area of skin to be treated
- ➔ For targeted treatment of the various sizes of skin, attach the relevant Aperture plate on the device surface.
- ➔ Attach Comb, when radiation exposure is meant for the patches over scalp or type of skin is sensitive.
- ➔ Carry out the therapy procedure on alternate day
- ➔ Wear the goggle to avoid direct exposure of NBUVB radiation to the eyes if exposure is not intended for this part.
- ➔ Close the eyes firmly, if intended area of treatment is exterior of eye lids or nearby area
- ➔ This is advisable to use dark room for effective exposure of of NBUVB radiation

***FOR Proper and effective treatment, refer the hand book thoroughly**

GENERAL PROTOCOLS FOR TREATMENT OF VITILIGO/PSORIASIS & OTHER SIMILAR CONDITIONS WHILE USING NBUVB DEVICE:

1. DOSE PROTOCOLS WITH RESPECT TO VARIOUS SKIN TYPE FOR THE TREATMENT OF VITILIGO:

	SKIN TYPE	Initial dose of NBUVB (Seconds)	Incremental Dose of NBUVB after each session (Seconds)	Max. dose of NBUVB allowed (Seconds)
	TYPE-1 VERY FAIR / PALE	40-Sec	10-Sec (40 + 50 + 60 + 70 + 80 90 + 100 + 110 + 120 +up to Erythema)	600-Sec
	TYPE-2 FAIR	60-Sec	15-Sec (60 + 75 + 90 + 105 + 120 + 135+ 150 + 165 +up to Erythema)	600-Sec
	TYPE-3 LIGHT BROWN / OLIVE	80-Sec	20-Sec (80 + 100 + 120 + 140 + 160 + 180 +up to Erythema)	600-Sec
	TYPE-4 OLIVE TO BROWN	100-Sec	25-Sec (125 + 150 + 175 + 200 + 225+ 250 +up to Erythema)	900-Sec
	TYPE-5 DARK BROWN	120-Sec	30-Sec (120 + 150+ 180 + 210 + 240 + 270 +up to Erythema)	900-Sec
	TYPE-6 DEEP PIGMENTED DARK BROWN	140-Sec	35-Sec (140 + 175 + 210 + 245 + 270 + 305 +up to Erythema)	900-Sec

- * Above treatment protocols are applicable to NBUVB device without attaching Comb
- * Comb should be attached in case of patches on scalp or in case of sensitive skin area / sensitive skin type
- * Dose at which Light pink skin colour is perceived, is known as **Erythema**. This erythemic condition defines the compatible dose for patient and that fixed compatible dose must be continued.

2. FREQUENCY OF NBUVB PHOTOTHERAPY USAGE:

Three times per week or alternate day

3. IF SUBSEQUENT DOSES ARE MISSED:

MISSED DURATION	CHANGES MADE IN PHOTOTHERAPY TIME
4 - 7 Days	Keep dose same
1 - 2 Weeks	Restart at last exposed dose - 25%
2 - 3 Weeks	Restart at last exposed dose - 50%
More Than 3 Weeks	Restart again right from scratch

4. ADAPTATION INCREMENTAL OF DOSE:

After prolonged compatible dose, due to adaptation principle, erythemic condition may be disappeared. In that case, again we should increase the compatible dose by 20%, until new erythemic condition is achieved.

5. PROBABLE COURSE OF NARROW BAND UVB:

- ➔ 18-36 session needed prior to assess treatment response
- ➔ ≥ 48 sessions needed before before discontinuing phototherapy due to lack of response
- ➔ ≥ 72 sessions needed before before stopping phototherapy especially in slow responders

6. DOSE ADJUSTMENT AFTER REPLACEMENT OF NEW LAMP:

Reduce dose by 20% and again verify the perceptible erythemic condition

7. COMPLETE EVALUATION OF THE TREATMENT:

Ideally 3 months can be considered for complete evaluation of the treatment along with your skin consultant

8. POST NB-UVB EXPOSURE:

Avoid additional sun exposure, or it is recommended to apply broad-spectrum sun screen regardless of skin type , minimum skin protection factor of (SPF) ≥ 30

9. RECOMMENDED EMOLIENT/TOPICAL OILS BEFORE NB-UVB SESSION IN PSORIATIC SKIN:

Following topical oils can be applied over the skin to minimize dryness and to enhance penetration of exposed NBUVB energy, such as elbows and knees, eventually enhances the efficacy of treatment and accelerate the clearance of Psoriasis. We can use either of the following options.

Mineral oil:

This is a hydrocarbon compound also known as paraffin, white mineral oil & Nujol

Baby oil:

This is one of the mineral oil, only difference is fragrance is added.

➔ Glycerin, Olive oil, Vaseline oil, Oleic acid 5% are also highly recommended

* **Note:** Coconut oil, sun flower oil are not recommended

10. SPECIAL SITES PRECAUTION DURING NBUVB SESSION:**During Face exposure:**

Cover the face during phototherapy, if not involved for treatment

Male Genitalia:

Use shields or cover the same by any means if not involved for treatment

Female Areola:

Use sun screen, shields or cover the same by any means

11. POST THERAPY SKIN CARE:

- ➔ Don't take immediate bath after NBUVB therapy
- ➔ Bath with mild soap/baby soap or cleanser
- ➔ Shorten shower times
- ➔ Limit to one shower per day
- ➔ Use warm (not hot) water only
- ➔ Moisturize skin twice daily (not immediate before NBUVB Phototherapy, may reduce efficacy)
- ➔ Best time to moisturize skin is after bathing

12. POST MAINTENANCE OF TREATMENT REGIM:

After recovery, it is not advisable to abandoned therapy on of sudden, rather should be minimize the days slowly.


- ➔ Dose, twice a week for initial 4-weeks
- ➔ Dose a week for next 4-weeks

that means total number of doses will be 12 , for post maintenance of treatment regim.

13. WHAT SHOULD I DO, IF MY SKIN BECOMES RED ?

Some reddening of the skin after photo therapy is normal. It will take between four to 24 hours after your photo therapy session for the skin to become red. However, if the skin looks very red, or is painful then you need to change your treatment schedule a little.

If your skin is red after your last light therapy session, apply moisturiser to your skin and use this diagram to decide what to do next.



14. COMBINATION THERAPY:

IMMUNOMODULATORS (IMMUNE SYSTEM SUPPRESSOR) + NBUVB THERAPY:


This is practically proven that in combination therapy, excellent results have been achieved in treatment of Vitiligo and Psoriasis. For better understanding, we need to know probable causes of Vitiligo and Psoriasis ailments. We will not go into anatomical detailed explanation rather go through simplified explanation and could be given as follow:

HOW VITILIGO IS TRIGGER:


**AUTO IMMUNE SYSTEM
IS DISTURBED**



**BIOCHEMICAL
CHANGES IN BODY**



**HEREDITARY OR
DISTURBED DNA
PATTERN**





**ENVIRONMENTAL
EFFECT**



**VIRUS, BACTERIA
OR PARASITES**



NORMAL SKIN:**Normal Pigmented skin****DEPIGMENTED SKIN / VITILIGO:****White Patches**

PATHOGENESIS OF VITILIGO:

Due to autoimmune disorder, T-Cells become overactive and attack the melanocytes. Melanocytes are the type of the skin cells, which are prime responsible to define the skin colour of epidermis (upper layer of skin). Due to absence of melanocytes, upper layer looks white, this condition is known as Vitiligo.

TREATMENT OF VITILIGO WITH COMBINATION THERAPY: IMMUNOMODULATOR, TOPICAL CREAMS/OINTMENTS:

Progressive nature / spreading of Vitiligo can be prevented by Immunomodulators. By using certain topical ointments or creams, T-Cells overactive nature can be controlled or minimize & eventually melanocytes can be retained.

Topical ointment, i.e: Tacrolimus (0.03% for face application) and Tacrolimus (0.1% for other body parts except genital parts) can be applied daily at night (because sun exposure should be avoided).

- At least 2-hrs gap (before and after therapy session) between NBUVB Phototherapy session and Topical application is recommended.

- This can be applied up to 16-weeks without any significant or serious side effects and also this does not include any type of steroids.

Therefore Tacrolimus ointment is known as one of the safest topical ointment.


NBUVB Therapy:

When NBUVB Phototherapy regime is followed as per earlier discussed exposure protocols, Vitiligo can be treated with good efficacy.

When desired skin is exposed with NBUVB radiation, melanocytes are stimulated, further melano blasts are proliferated and migrated, which results in repigmentation of skin.

Hence combination therapy, Immunomodulators (Immune system suppressor) + NBUVB Therapy brings synergetic action in the treatment of Vitiligo / Leucoderma. Eventually fastest recovery results are achieved

PATHOGENESIS OF PSORIASIS:



PSORIASIS:

Due to malfunction of autoimmune system, T-Cells react differently in such a manner that the skin cell are produced rapidly. Precisely speaking, in healthy skin, normal cycle of skin cell production is 28-30days, but in Psoriatic condition, skin cell production is 6-8 days. Thus don't allow enough time for the older cells to be dried up and new cells get piled up, eventually upper layer is appeared to be swollen, scally with silvery flakes and red blister.

TREATMENT OF PSORIASIS WITH COMBINATION THERAPY: IMMUNOMODULATORS (IMMUNE SYSTEM SUPPRESSOR) + NBUVB THERAPY:

Spreading of disease can be prevented by immunomodulator. By using certain topical ointments or creams, T-cell overactive nature can be controlled or minimize and eventually abnormal cell producing cycle can be controlled.

Topical ointment, i.e: Tacrolimus (0.03% for face application) and Tacrolimus (0.1% for other body parts except genital parts) can be applied daily at night (because sun exposure should be avoided).

At least 2-hrs gape (before and after therapy session) between NBUVB Phototherapy and Topical application is recommended. This can be applied up to 16-weeks without any significant or serious side effects and also this does not include any type of steroids. Therefore Tacrolimus ointment is known as one of the safest topical ointment.

When NBUVB Phototherapy regime is followed as per earlier discussed exposure protocols, Psoriasis can be treated with good efficacy. When desired skin is exposed with NBUVB Phototherapy device, T-Cells are deactivated, which results in prevention and control of premature production of unnatural cycle of skin cells.

Hence combination therapy, Immunomodulators (Immune system suppressor) + NBUVB Therapy brings synergetic action in the treatment of Psoriasis and ailment could be managed very well.

* Note: Though there is availability of so many types of the topical creams/ointments in the market, can be used as combination therapy, but highly recommended to consult your skin consultant or relevant Medical professional.

FEW GOLDEN TIPS:

HOW TO USE NBUVB PHOTOTHERAPY DEVICE?

- ➔ During NBUVB Phototherapy session, Device should be kept as close as possible against white patch, even you may touch the device on white patch & must be carry out thrice a week or alternate day.
- ➔ As a thumb rule (in case of any type of skin), you should start the session for 60 Seconds on 1st day and during next session phototherapy should be carried out for 75 Seconds.
- ➔ During each alternate session, go on increasing the exposure timing @ 15 seconds (60 sec + 75 sec + 90 sec + 105 sec + -----) until white patch seems to be light pink colour, this is called erathemic condition.
- ➔ For example light pink colour is perceived at 120 seconds, then that particular time will be compatible dose for your skin, will result in fast recovery from white patch
- ➔ Now you should expose the patch for 120 seconds for forthcoming each and every session, eventually will lead to fast recovery.

FEW TIPS TO GET GREAT RESULTS:

- ➔ Apply suitable moisturizer to skin after phototherapy to retain good tone and texture of the skin
- ➔ Avoid excess water during bath
- ➔ Try to use luke warm water

- ➔ Try to use baby soap or mild body wash liquid
- ➔ Avoid direct sun light during out door activities or you may use sun screen having at least SPF:30
- ➔ Take chewable tablet of calciquick D3, 60000 IU, by Dr Morepen, Dosage: once a week, can be taken regularly (once a week basis) up to 6-months. Discontinue for month and again continue and so on,this will expedite the recovery rate @21% from Vitiligo
- ➔ This is recommended to take NBUVB Phototherapy session in dark room, otherwise spurious / ambient light may interfere to to desired radiation of exposed energy
- ➔ NBUVB Phototherapy device, 311nm should be used regularly with accurate timing.
- ➔ Always follow the instructions of your Doctor / Dermatologist, because he knows you and skin condition up to the mark
- ➔ For fast recovery and get rid of ailments, always make strong positive affirmations to your subconscious mind. Really it works!


Disclaimer:

This hand book is not designed to enable any one to undertake any kind of phototherapy treatment independently, rather highly recommended to get appropriate supervision by suitably qualified medical professional. Guidelines and & treatment plan is suggested plan only. Medical professionals are free to adapt, amend and reuse the hand book contents at their own way.

Neither the manufacturer nor the trial based guidelines described in handbook will accept any responsibility or liability for any harm or side effect that may arise from the use of this device.

*** SPECIFICATIONS CAN BE CHANGED WITHOUT PRIOR NOTIFICATION DUE TO CONTINUOUS R&D.**

Manual/Hand Book




NBUVB, 311nm Lamp 5-Star Rated CE Certified/PHILIPS (Made in Poland)



NBUVB Device




NBUVB Device Container



UV Protective Goggle with Case & Micro Fiber



PHOTOMETRIC RESPONSE OF NBUVB DEVICE:



Eczema



Psoriasis



Vitiligo



DEALERS DETAILS :

Warranty: 1-Year
(Except Lamp)



MFG. BY:


VISION TECH[®]
ELECTRONICS

C/9, Bansidhar Estate
Opp. Suprabhat Industrial Estate
O/s Dariyapur Gate, Bardolpura
Ahmedabad (Gujarat, India) - 380004
M:9925799688



E-mail: info@visiontechproducts.com

Web Site: www.visiontechproducts.com

 **YOUTUBE** <http://www.youtube.com/@visiontechproducts>